

---

# HowToPulsePicker Documentation

*Release 1.0.0*

**[gabriele.giovanetti@xfel.eu](mailto:gabriele.giovanetti@xfel.eu)**

**Dec 16, 2021**



---

## Contents

---

<b>1</b>	<b>How to run scans with the Pulse Picker</b>	<b>3</b>
1.1	Pulse picker specific devices . . . . .	3
1.2	Other devices needed for Pulse Picker scans . . . . .	4
<b>2</b>	<b>Example setup for Pulse Picker scans</b>	<b>5</b>
2.1	Devices configuration remarks . . . . .	5
<b>3</b>	<b>Indices and tables</b>	<b>7</b>



Contents:



---

## How to run scans with the Pulse Picker

---

### 1.1 Pulse picker specific devices

#### 1.1.1 The Pulse Picker

The Pulse Picker is a hardware device that controls a shutter. It is designed to let a single pulse out of a train pass through (no matter which pulse). It's currently installed in SPB. From the Karabo point of view it looks like a bunch of Beckhoff digital inputs/outputs. We use two of them to control it: We set SPB\_XTD9\_PPU/DCTRL/TRIG *ON* to start the pulse picker (and it goes back to *OFF* when it is done) We read back SPB\_XTD9\_PPU/DCTRL/STEP2 as its raising edge is the best approximation we have of the time of the passing pulse. *Note:* motor SPB\_XTD9\_PPU/MOTOR/ROTOR is used behind the scenes in follow-target mode. Do not mess with its settings. Just in case, it may help to know that slit is open when motor actual position is  $(2 * N) * 30$  for any N, closed when it is  $(2 * N + 1) * 30$ . So resting value for actual position should be any odd multiple of 30, and *stepSize* value should be 30.

#### 1.1.2 PulsePickerTrigger

It is a middlelayer device. It repeatedly activates Pulse Picker for a given *Acquisition Time* at a given *Repetition Rate* (thus letting *Number of Pulses* pass through). As a pulse passes through, the device writes on and output channel a boolean with the same train Id of the pulse (as described in *Pulse Picker*). It is meant to be *Karabacon*-friendly.

When the “Pass-Through” option is enabled, it opens the pulse picker slit (assuming it is closed at the beginning), keep it open for at least *acquisitionTime*, then closes it again. At the same time *ImagePicker* is disabled.

In both cases if averagers names are set, it takes care of resetting averaged scan values.

#### 1.1.3 ImagePicker

Yet another device from the *imageProcessor* family. It gets an image stream in input (e.g. from a camera), and the train id on a second one (as provided from *PulsePickerTrigger*). It writes on an output channel only the images whose train id matches the one received by *PulsePickerTrigger*. Users can set an offset value to delay or anticipate the image

train id with respect to the pulse one. The ImagePicker can be disabled by setting to `True` the **Disabled** boolean. This way the images in input are simply copied to output, as if the ImagePicker was not present.

## 1.2 Other devices needed for Pulse Picker scans

### 1.2.1 LimaBaslerCamera

Used to acquire beam image. It is supposed to acquire at 10 Hz.

### 1.2.2 ImageProcessor

Used to evaluate beam width ( $4 * \sigma_X$ ) and height ( $4 * \sigma_Y$ ) on the relevant images (*ImagePicker* Output).

### 1.2.3 GenericAverage

This device is used to average beam width and height read from *ImageProcessor* over values.  $N$  must be  $\leq$  *PulsePickerTrigger* *Number of Pulses*. In the ideal case (i.e. all pulses lead to a valid value) they are equal. Devices and properties to be averaged can be defined at configuration time. Output averaged values have the form of nodes with *average* and *nOfValues* properties. The name of the node can be configured, **Note that such names must match the hard-coded interface of Karabacon** Average values have the form of nodes, whose name can be selected at con

### 1.2.4 Karabacon

Meant to scan over a motor position, use *PulsePickerTrigger* as trigger and *GenericAverage* as a data source.

---

## Example setup for Pulse Picker scans

---

The following informations are updated to 11-10-2018. They can anyway be useful for refernece as an example.

In order to test pulse picker scans the following system was set up. It may be found in subproject *SCANTOOL* in domain *SPB*.

Device Class	Device ID
LimaBaslerCamera	MY_FAVOURITE_CAMERA
ImagePicker	SPB_XTD9_PPU/PROC/IMPICK
ImageProcessor	SPB_XTD9_PPU/PROC/IMPROC
PulsePickerTrigger	SPB_XTD9_PPU/MDL/TRIGGER
GenericAverage	SPB_XTD9_PPU/MDL/BW_AVG
Karabacon	SPB_KARABACON

### 2.1 Devices configuration remarks

- ImagePicker SPB\_XTD9\_PPU/PROC/IMPICK
  - **Input image:** connected output channel = MY\_FAVOURITE\_CAMERA:output  
*The camera name here is all you need to change when switching image soource*
  - **Input train id:** connected output channel = SPB\_XTD9\_PPU/MDL/TRIGGER:output
- ImageProcessor SPB\_XTD9\_PPU/PROC/IMPROC
  - **Input:** connected output channel = SPB\_XTD9\_PPU/PROC/IMPICK:output
- PulsePickerTrigger SPB\_XTD9\_PPU/MDL/TRIGGER
  - **averageNames** = [ 'SPB\_XTD9\_PPU/MDL/BW\_AVG' ]
  - **repetition rate** and **acquisition time** can be configured to achieve the desired number of pulses per scan point
- GenericAverage SPB\_XTD9\_PPU/MDL/BW\_AVG

– **properties table**

deviceName	propertyName	averageOutputName
SPB_XTD9_PPU/PROC/IMPROC	beamWidthId	foursigma <sub>x</sub>
SPB_XTD9_PPU/PROC/IMPROC	beamHeightId	foursigma <sub>y</sub>

**Note:** averageOutputName values can be changed at configuration time **but** they have to match the hardcoded values in Karabacon interface, so do not change them if you want to perform a scan.

- **Window Size:** it could make sense to have this matching the PulsePickerTrigger numberOfPulses but it is not strict

- Karabacon SPB\_KARABACON

- Main control scene:

- \* Motor environment: the motor you want to scan on (axis: default)
- \* Data sources:

Alias	DeviceId	value
[e.g.] beamwidthX	SPB_XTD9_PPU/MDL/BW_AVG	foursigma <sub>x</sub>
[e.g.] beamwidthY	SPB_XTD9_PPU/MDL/BW_AVG	foursigma <sub>y</sub>

- \* Trigger Sources: put here the device id of PulsePickerTrigger (i.e. SPB\_XTD9\_PPU/MDL/TRIGGER)

## CHAPTER 3

---

### Indices and tables

---

- `genindex`
- `modindex`
- `search`